

THE SAVI® BREAST BRACHY THERAPY APPLICATOR: Partial Breast Radiation Therapy for Early-Stage Breast Cancer



WHAT IS SAVI?



The Strut Adjusted Volume Implant (SAVI) applicator delivers a form of accelerated partial breast radiation (APBI) known as breast brachytherapy, which delivers radiation from inside the breast. This allows physicians to precisely deliver treatment to the tumor cavity and surrounding tissue, and offers two key advantages: it spares healthy tissue from unnecessary radiation and reduces treatment time from 6 weeks to just 2–5 days*.

HOW DOES IT WORK?


SAVI uses a bundle of catheters to deliver radiation. The catheters are expanded to fit the size and shape of the cavity, and a tiny radioactive seed travels through each catheter. This allows every catheter to deliver an individual dose of radiation, so physicians can customize treatment for every patient, regardless of the size, shape, or location of their tumor cavity.

BENEFITS OF SAVI BREAST BRACHY THERAPY


- Customizes radiation to individual patient
- Excellent local control of cancer
- Reduces toxicity from unnecessary radiation
- Treatment lasts just 2–5 days*

*based on clinicians' preferred treatment protocol. Khan, Atif J. et al. Three-Fraction Accelerated Partial Breast Irradiation (Apbi) Delivered With Brachytherapy Applicators Is Feasible And Safe: First Results From The Triumph-T Trial. International Journal of Radiation Oncology. 2019. Online publication.


1 PLACEMENT
The Physician gently inserts SAVI in a closed position through a small incision, and then expands the catheters to conform to the shape of the cavity. The ends of the catheters remain outside the breast for radiation delivery. SAVI remains in the breast during the entire course of treatment.



2 CUSTOMIZED TREATMENT PLANNING
The radiation oncologist performs a CT scan and uses those images to create a customized treatment plan to deliver radiation to the areas that need it most, while protecting the skin, chest wall, heart and lungs.



3 RADIATION DELIVERY
Treatment is delivered twice a day for 2–5 days: Each treatment lasts approximately 10 minutes. The ends of the catheters are connected to a device—called an HDR afterloader—and a computerized system painlessly delivers a wire with a tiny radioactive seed into each catheter.



The seeds are completely removed after each treatment—no radiation remains in the patient's body between treatments. After the last treatment, the physician closes the device and removes it through the same incision in which it was inserted.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



merit.com
ciannamedical.com

©2019 Cianna Medical, Inc. All rights reserved. Patents and trademarks listed at www.ciannamedical.com/patents

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733